

INTRODUCTION

The complete blood count (CBC) is the most commonly ordered blood test in a healthcare set up. The fully automated sophisticated cell counters give us very accurate and precise blood counts rapidly along with a number of additional newer parameters which give valuable extra information added to the routine CBC. The origin of the modern cell counters date back to 1953, when Wallace H Coulter, an American engineer by training, obtained a patent for a principle for counting and sizing microparticles suspended in a fluid.¹ The principle was named the 'coulter principle' and is the basis of the technology of most cell counters even today.

A basic CBC comprises of the following parameters

Haemoglobin (Hb)

Haematocrit (Hct)

Red blood cells(RBCs) and indices

Total Red Cell count (TRBC)

Mean corpuscular volume (MCV)

Mean corpuscular haemoglobin (MCH)

Mean corpuscular haemoglobin concentration (MCHC)

Red cell distribution width (RDW)

White blood cells (WBCs)

Total white cell count (WBC)

White cell count differential (DC)- 3 parts /complete 5 part differential count.

Absolute white cell counts

Platelets

Platelet count (PC)

Mean platelet volume (MPV)

Platelet distribution width (PDW)

Over and above all these, few new parameters available on some analyzers have found clinical utility and are being incorporated into the routine CBC. Two new US FDA approved parameters are the Reticulocyte Haemoglobin Content (CHr) or Reticulocyte Haemoglobin Equivalent (Ret-He) and Immature Platelet Fraction (IPF). Promising research parameters include Leucocyte Positional Parameters(LPP).²


The Hb and the Hct

The haemoglobin is measured after the RBCs are lysed and the Hb converted to a stable coloured compound which is measured spectrophotometrically. Hence, like lipemia and hyperbilirubinemia can interfere with the results.

The Hct can be used synonymously with the Hb to classify a sample as anaemic or polycythemic. The Hct is the given volume of packed RBCs as compared to the plasma and can be measured by a microhaematocrit centrifuge (as in Sysmex analyzers) or calculated (as in Beckman Coulter instruments) by taking into account the Hb and the MCV. Due to differences in the methodologies there can be considerable variation in the Hct between analyzers while the Hb value is more reproducible. Hence, it is better to use the Hb to make clinical decisions.

TRBC and MCV

The most common technology for counting and sizing the blood cells is the impedance technology based on the coulter principle. The basic counting chamber comprises of two electrodes with a partition between them. The partition has an aperture through which electrolyte flows and maintains flow of electrical current between the electrodes (Figure 1). The cells are suspended in the electrolyte and when an individual cell passes through the aperture, there is a break in the flow of current as blood cells are poor conductors of electricity. Each break is counted as a pulse. The number of pulses generated give the cell count. The TRBC and the platelet counts are obtained in this manner.



3-part Diff technology

V0706

Fig. 1: Impedance based cell counter.
Source-<https://www.sysmex.com>




Fig. 2: CBC printout




Fig. 3: RBC Histogram

The height/size of a pulse generated depends on the size of the cell causing it. In the case of RBCs, the average height of the pulses gives us the MCV. The MCV is expressed in femtolitres (fl) and red cells with MCVs in the range of 80-100fl are considered normocytic. Cells with MCV less than 80fl are called microcytic. An MCV above 100fl would classify the cells as macrocytic. Besides numerical values, the instruments also give RBC volume distribution histograms.

MCH AND MCHC

Having obtained the Hb, Hct, TRBC and MCV, the instruments give us the calculated parameters of MCH ($Hb \times 1/TRBC$) and MCHC ($Hb \times 1/Hct$) which give information about the concentration of Hb in the RBCs. Values below the reference range are used to classify the RBCs as hypochromic. Samples with higher concentrations of MCH indicate dense RBCs as may be seen in megaloblastic anaemias or in the presence of spherocytes. Very high (non-physiological) values of MCV, MCH and MCHC indicate presence of red cell clumping, often due to cold agglutinins. The abnormalities are reversed on reanalyzing the samples after incubation at 37 deg C. The CBC printout with the characteristic RBC indices is diagnostic of red cell agglutination. The following RBC histogram illustrates the same (Figure 2).

RDW

From the RBC volume distribution histogram, the analyzers calculate an index of variation of RBC size (anisocytosis) called the RDW. Anisocytosis is appreciated on the peripheral smear (PS) but it is not possible to




Fig. 4: Algorithmic approach to Anaemia

quantify anisocytosis and give a numerical value. On the other hand, on the analyzers, a numerical value of the RDW is obtained as a percentage coefficient of variation (CV%) or as the standard deviation of the mean size (SD). A raised RDW indicates a heterogeneous population of RBCs from the size point of view.³ In nutritional anaemias, this parameter is often the first to be altered and a raised RDW even with a normal MCV warrants a PS examination to look for a small population of microcytic or macrocytic RBCs. It is important to realise that the MCV is just an average value and even if a population of microcytes /macrocytes are present in a sample, the MCV may fall in the normal range if majority of the cells are normocytic. Moreover, an admixture of microcytes and macrocytes as may be seen in a case of nutritional anaemia may average out the MCV, causing it to fall in the normal range. However, the RDW would be much raised. Thus, the RDW is an useful parameter to differentiate between nutritional and other causes of anaemias.

The following RBC histogram shows marked anisocytosis and a raised RDW (Figure 3).

RETICULOCYTE COUNTS

An underutilised but useful parameter in the work up of anaemias is the reticulocyte count. As reticulocytes are the youngest forms of RBCs, their presence in the blood establishes the fact that the bone marrow is producing RBCs. Increased bone marrow production of RBCs is reflected in a raised reticulocyte count. Unlike mature RBCs, the reticulocytes contain remnant RNA which is stained by a supravital dye and which appears as a 'blue reticulum' under the microscope.⁴ Thousand RBCs have to be counted and reticulocytes enumerated and expressed as a percentage of RBCs which can be difficult and unreliable. Today's cell counters provide automated




Fig. 5: Scattergram from a Sysmex X class analyzer. Source- <https://www.sysmex.com>

reticulocyte counts on the basis of their larger volumes and presence of RNA.

For ease and familiarity, the reticulocytes continue to be expressed as a percentage of the RBCs which may lead to falsely raised values when the absolute number of RBCs is low as seen following haemolysis or blood loss. A better way is to 'correct' the reticulocyte count for the patient's Hb or Hct. A simple formula is- Corrected reticulocyte count = % of reticulocyte \times patient's Hct/45.⁵

Getting used to absolute reticulocyte numbers provided by the automated analyzers may be a better solution but limited by the limited access to automated reticulocyte counts.

CBC IN THE WORK UP OF ANAEMIAS (FIGURE 4)

The following is a simple algorithm for the start up of the work up of a case of anaemia using the Hb, red cell indices and reticulocyte count.⁵

If the reticulocyte count is raised and blood loss is ruled out, few routine biochemical tests like serum bilirubin and fractions, haptoglobin, lactate dehydrogenase (LDH), etc, may establish ongoing haemolysis. Added to these, a peripheral smear examination may help to clinch the exact diagnosis/cause of haemolysis. For example, finding a spherocyte would suggest an autoimmune haemolytic anaemia or hereditary spherocytosis while a 'bite cell' or 'blister cell' would suggest oxidant damage as seen in G6PD enzyme deficiency.

Alternatively, one may look at the RBC indices. Macrocytic anaemias may be due to B₁₂/folate deficiencies but may be also seen with use of anti retroviral drugs, liver disease, hypothyroidism.⁶ Macrocytic anaemia as a part of generalised pancytopenia may be seen in megaloblastic anaemia or aplastic anaemia. A unique type of myelodysplastic syndrome, associated with loss of the long arm of chromosome 5 (5q- syndrome) is

characterised by macrocytic anaemia and thrombocytosis. Increased reticulocytes seen in haemolytic anaemias may lead to slight macrocytosis. Normocytic indices may be seen in ACD, haemolytic anaemias and even combined nutritional anaemias.

In the presence of microcytosis, the TRBC (raised in thalassemia) and RDW (increased in IDA) can help differentiate between iron deficiency and beta thalassemia trait. Indices like the Mentzer's index⁷ have been put forward to make this differentiation more objective. If the quotient of the MCV (in fl) divided by the TRBC (in millions/cubic mm) is less than 13, thalassemia is more likely. If the result is greater than 13, iron deficiency is more likely. Co-existing iron deficiency and beta thalassemia trait make these values redundant. Beta thalassemia carriers have microcytic hypochromic RBC indices with a normal RDW as the red cells are uniformly small, i.e., homogeneously microcytic. On the other hand, iron deficiency shows microcytosis along with an altered RDW.

CHr/ Ret-He

A close differential diagnosis of iron deficiency anaemia (IDA) is anaemia of chronic disease (ACD), characterised by functional iron deficiency (FID) where iron is not incorporated into developing erythroid cells in spite of adequate bone marrow iron stores. It is now well established that release of hepcidin from liver cells under the influence of inflammatory cytokines like interleukin -6 is responsible for blockage iron iron movement from iron storing macrophages into the developing red cells.⁸ This block also occurs at the level of the intestinal cells wherein absorbed iron cannot enter the blood stream. Over weeks, this lead to formation of RBCs with low MCV and low MCH. The traditional biochemical tests like serum iron, transferrin and ferritin levels fail to convincingly distinguish between IDA and ACD. An improvement over these tests is the soluble transferrin receptor level (sTfr) which is raised only in IDA and remains low in ACD.⁹ This parameter is limited by its complexity of testing and its relative unavailability.

Changes in the RBC parameters as described occur over the life span of RBCs (120 days) while changes occurring in reticulocytes occur quickly, i.e., within 2-3 days (the life span of the reticulocytes). CHr provided by instruments manufactured by Siemens and Ret-He provided by Sysmex analyzers measure reticulocyte haemoglobin content or its equivalent directly. Both are helpful parameters to detect iron insufficient erythropoiesis,¹⁰ either due to pure and simple iron deficiency or due to FID at a very early phase. Effective treatment of these conditions bring about a rapid change or increase in the CHr or Ret-He as these are measurements on the reticulocytes. Moreover, this information is obtained at the cost and at turn around time of a routine CBC. However, it should be noted that the Ret-He is a size based assay and will be useless in the presence of an underlying beta thalassemia or in the presence of concomitant megaloblastic anaemia.

It has been shown that using a combination of sTfr, serum




Fig. 6: Platelet clumps in EDTA

ferritin and the Ret-He, it is possible to accurately diagnose and provide guidelines for appropriate management of all microcytic hypochromic anaemias (Thomas plot, 11). Both Ret-He and CHR have been incorporated into guidelines for management of anaemia in chronic renal disease patients and FID.¹²

A high haematocrit which is inexplicable (no history of living at high altitude, no smoking, no chronic obstructive airway disease, etc) should prompt one to work up for a primary haematological disorder like polycythaemia vera (PV). A the Jak2 mutation performed on peripheral blood is warranted.¹³

WBCS

Popular technologies for enumeration and differentiation of WBCs include the Volume Conductivity and Scatter (VCS) used by the Beckman Coulter instruments (from Beckman Coulter Inc, USA) and the Fluorescence Flowcytometry (FFC) principle used by the Sysmex analyzers (from Sysmex Corporation, Japan).

The VCS principle differentiates cells based on their volume measured by impedance, conductivity of radio frequencies which give information about the internal structure of cells and and scatter of light which tells us about surface structure and granularity. In the FFC technology, cells are hit by a laser beam which gets scattered in various directions. The forward scattered light is a measure of the cell volume while the side scattered light is a measure of cell granularity. When a fluorescent dye is used, side fluorescence gives information about the nucleic acids in the cells. Thus, monocytes have the highest forward scatter owing to their large size while eosinophils have the highest side scatter owing to the presence of the prominent granules. immature granulocytes have high fluorescence. Thus, various scatter plots are generated and the total WBC and the differential counts obtained.

Leucocyte positional parameters- Numerical data coordinates generated by VCS technology are available as research parameters on LH Beckman Coulter analysers. Volume, conductivity and light scatter for each WBC type

are characteristic. Any deviation from these as given by the mean value and standard deviation of each of these parameters reflect changes in size and complexity for a particular cell type which in turn may be reflective of specific diseases. These positional parameters have been used to differentiate between chronic lymphoproliferative disorders (CLPD) like chronic lymphocytic leukaemia (CLL) vis a vis reactive lymphocytosis as in the former, the cells are uniformly small.¹⁴ Changes in the parameters of neutrophils have been found to indicate sepsis and bacterial infections.¹⁵ Changes involving lymphocytes and monocytes have also been used to flag for presence of malaria.¹⁶ As of now, these are still research parameters and yet to be incorporated into routine clinical practice.

The CBC will indicate if there is leucocytosis or leucopenia. The type of cell involved suggests the possible differential diagnoses. For example, increase in neutrophils could suggest a bacterial infection, an inflammatory state, steroids use, etc. An extremely high count with left shift suggests a chronic myeloid leukaemia (CML) or a leukamoid reaction. A bimodal peak of myelocytes and neutrophils alongwith presence of basophils point towards CML and the molecular test for the diagnostic fusion gene, BCR:ABL1 should be asked for.

Similarly, finding sustained absolute lymphocytosis without any chronic infection would warrant a haematological work up for a CLPD like CLL. Lymphoma spillover into blood may also cause a lymphocytosis. Peripheral smear may reveal characteristic morphology but immunophenotyping using flowcytometry is the main stay of diagnosis.¹⁷

Eosinophilia is most often secondary to allergic conditions, parasitic infestations, drugs, vasculitides. If no such cause is found after reasonable investigations and if eosinophilia is sustained, a primary eosinophilia (clonal or hypereosinophilic syndrome) needs to be considered. Molecular studies for FIP1L1-PDGFR mutation are important to determine response to a tyrosine kinase inhibitor like Imatinib.¹⁸

Monocytosis may be seen in some chronic infections or during recovery of counts following chemotherapy or drugs. However, sustained absolute monocytosis of at least 1000/ul in absence of chronic infections warrant investigating for a haematological disorder like chronic myelomonocytic leukaemia in an adult and juvenile myelomonocytic leukaemia in children.

In the presence of leucopenia, it is prudent to look at the absolute cell counts to see if any particular cell line is predominantly suppressed. Neutropenia is most commonly encountered due to drugs and severe infections which would be apparent from the patient's clinical condition. When neutropenia occurs alongwith anaemia and thrombocytopenia, a serious bone marrow disorder like aplastic anaemia should be thought of. Inherited conditions of neutropenia are rare but known (e.g., Kostman syndrome). Acquired lymphopenia may be seen following use of monoclonal anti CD20 antibodies




Fig. 7: RBC & Platelet Histogram

like Rituximab.¹⁹ Absolute monocytopenia is a feature of Hairy cell leukaemia while eosinopenia is a feature of sepsis.

PLATELETS

A small percentage of people develop platelet clumping when blood is collected in EDTA (the preferred anticoagulant for CBC analysis) because they carry antibodies against certain epitopes of platelets unmasked in the presence of EDTA (Figure 6).²⁰ In such cases, a repeat count after a repeat collection in citrate is indicated.

It should be noted that in the presence of platelet clumps it is not possible to give an accurate platelet count, whatever be the technology used.

The most common method of counting platelets is impedance as described earlier in the section on RBCs. The difference in size between the RBCs and platelets makes it possible to separate out the the platelets and RBCs. Problems occur when there are very large platelets which may be falsely picked up as RBCs leading to a falsely lowered platelet count.

Pathological conditions associated with large platelets like Bernard Soulier syndrome and familial platelet disorders with myosin heavy chain9 (MYH9) gene mutation are well documented but are rare.²¹ More commonly, in the eastern part of our country, a large chunk of the population have mild thrombocytopenia and very large platelets (constitutional macrothrombocytopenia) which cause further artefactual lowering of the platelet count on impedance cell counters.²² There is no risk of bleeding but the low platelet count leads to panic and unnecessary investigations.

Extremely microcytic or fragmented RBCs may be wrongly counted as platelets causing an inaccurate higher platelet count. This becomes more significant when the platelet counts are low.

In these situations, the platelet and the RBC histograms overlap with each other and the instrument gives a flag which warrants a peripheral smear check. A manual platelet count is useful in such situations.




Fig. 8: Measurement of IPF using platelet-F technology.
Source-<https://www.sysmex.com>

The following figure shows overlapping RBC and platelet histograms (Figure 7).

As manual platelet counts are difficult and time consuming, there have been attempts to innovate technology to overcome the problems with impedance counts like optical counts based on light scatter, use of monoclonal anti platelet antibodies (Abbott), fluorescent platelet counts (platelet-F, Sysmex),²³ etc.

Platelet -F involves use of a fluorescent dye which is highly specific for platelet membrane and counts obtained by this method have been found to be very accurate.²⁴

Besides the counts, platelet indices like the MPV and PDW have also been useful reflecting the average platelet size and platelet anisocytosis respectively.

A PS confirmation of thrombocytopenia is a must to rule out EDTA induced agglutination, presence of large platelets, fibrin strands, etc.

Infections and sepsis alongwith drugs are the commonest causes of thrombocytopenia encountered in a hospital set up wherein the clinical history is indicative of the cause. Thrombocytopenia associated with abnormalities in other cell lines may indicate a bone marrow disorder like an acute leukaemia. Isolated thrombocytopenia in the absence of any other cytopenia (except perhaps for an IDA secondary to bleeding) may be due to an immune thrombocytopenia (ITP) which may be primary or secondary to an associated immunological disorder.²⁵ ITP is a diagnosis of exclusion. It is crucial to rule out a thrombotic thrombocytopenia purpura (TTP) in the appropriate clinical setting by looking for fragmented cells (schistocytes) in the PS and analysing the serum LDH as TTP has a 90% mortality if plasmapheresis is not started early. Another life threatening situation is Heparin induced thrombocytopenia which needs to be considered in patients exposed to unfractionated heparin.²⁶

Thrombocytosis(>450,000/ul) is seen in inflammatory

332 states, infections and IDA. However, persistent thrombocytosis in the absence of these conditions may indicate an underlying chronic myeloproliferative neoplasm like CML, essential thrombocytosis, idiopathic myelofibrosis, etc. In such situations, bone marrow trephine biopsy and molecular tests for BCR;ABL1, JAK2 mutation, CALR mutations are warranted.²⁷

IMMATURE PLATELET FRACTION (IPF)

Akin to the reticulocyte count in RBCs is the IPF in the platelets (Figure 8). In the higher end Sysmex analyzers it is possible to reliably quantify the fraction of the most immature platelets containing RNA (reticulated platelets) which is expressed as the IPF % (reference range-1.1 to 6.1%).²⁸

In a case of thrombocytopenia, a high IPF indicates adequate bone marrow production with peripheral loss/destruction of platelets. It is seen that with recovering platelet counts, the IPF value falls. Thus, it is a very useful parameter in the diagnosis of various thrombocytopenias associated with increased peripheral destruction, eg., ITP/TTP.²⁸

IPF has been found useful in the management of Dengue. In a small study of 32 patients, 93.75% of the patients showed platelet recovery within 24–48 h if the IPF was more than 10% and 100% patients showed a recovery within 24 h of the fall of the IPF compared with the previous days.²⁹

SUMMARY

The CBC is often the starting point in the investigation of a sick patient and important clinical decisions are made based on the findings. Today, the CBC is obtained from automated analyzers which give fast, accurate and reproducible results. Parameters like the RDW are obtained only from the automated counters. Newer parameters like Ret-He have been incorporated into international guidelines for investigation of FID. A basic understanding of the technologies makes one aware of the pitfalls and possible problems. While interpreting an abnormal CBC, it is important to take the patient's clinical status into account and to review previous reports if any. For eg., neutrophilia in a hospitalized patient occurring over 2-3 days is likely to be due to an acute stress/infection while neutrophilia sustained over months in the absence of any apparent illness may be due to a myeloproliferative neoplasm like chronic neutrophilic leukaemia. Platelets are the most difficult cells to count due to the inherent problems associated with impedance technology. Methods like platelet-F give highly accurate counts. Newer parameters like IPF have been found to be in management of dengue and are being incorporated into the routine CBC.

REFERENCES

1. 2016 Wallace H. Coulter Foundation. <http://www.whcf.org>
2. Briggs C. Quality counts:new parameters in blood cell counting. *Int J Lab Hem* 2009; 31:277-297.
3. Lin CK, Lin JS, Chen SY et al. Comparison of haemoglobin

and red blood cell distribution width in the differential diagnosis of microcytic anaemias. *Arch Pathol Lab Med* 1992; 116:1030-32.

4. Bain B J, Lewis M S, Bates I 2006. Basic haematological techniques in Lewis, Bain and Bates (ed) Dacie and Lewis Practical Haematology. Elsevier Ltd.
5. Glader B. 2004. Anaemia:general consideration in Greer JP, Foerster J, Rodgers G M, et al (ed) Wintrobe's clinical haematology. Lippincott Williams and Wilkins.
6. Aslinia F, Mazza JJ, Yale SH. Megaloblastic anaemia and other causes of macrocytosis. *Clin Med Res* 2006; 4:236-241.
7. Mentzer WC. Differentiation of iron deficiency from thalassemia trait. *Lancet* 1973; 1:449-52.
8. Weiss G, Goodnough L T. Anaemia of chronic disease. *NEJM* 2005; 352:1011-1023.
9. Punnonen K, Irjala K, Rajamaki A. Serum transferrin receptor and its ratio to serum ferritin in the diagnosis of iron deficiency. *Blood* 1997; 3:1052-1057.
10. Brugnara C, Schiller B, Moran J. Reticulocyte haemoglobin equivalent (Ret He) and assessment of iron deficient states. *Clin Lab Haem* 2006; 28:303-308.
11. Thomas C and Thomas L. Biochemical markers and haematologic indices in the diagnosis of functional iron deficiency. *Clin Chem* 2002; 48:1066-1076.
12. Thomas Wayne D, Hinchliffe Rod F, Briggs C et al. Guideline for the diagnosis of functional iron deficiency. *British Journal of Haematology* 2013; 161:639-648.
13. Tefferi A, Gilliland DG. The JAK 2V617F tyrosine kinase mutation in myeloproliferative disorders:status report and immediate implications for disease classification and diagnosis. *Mayo Clin Proc* 2005; 80:947-958.
14. Silva M, Fourcade C, Fartoukh C et al. Lymphocyte volume and conductivity indices of the haematology analyzers Coulter GEN.STM in lymphoproliferative disorders and viral diseases. *Clin Lab Hem* 2006; 28:1-28.
15. Chaves F, Tierno B and Xu D. Quantitative determination of neutrophil VCS parameters by the Coulter automated hematology analyzer:new and reliable indicators for acute bacterial infection,. *Am J Clin Path* 2005; 124:440-444.
16. Briggs C, Da Costa A, freeman L, et al. Development of an automated malaria discriminant factor using VCS technology. *Am J Clin Path* 2006; 126:691-698.
17. Hallel M, Cheson BD, Catovsky D, et al. Guidelines for the diagnosis and treatment of chronic lymphocytic leukaemia. *Blood* 2008; 111:5446-5456.
18. Pardanani A, Brockman SR, Paternoster SF et al. FIP1L1-PDGFR α fusion:prevalence and clinicopathologic correlates in 89 consecutive patients with moderate to severe eosinophilia. *Blood* 2004; 104:3038-3045.
19. Plosker GL, Figgitt DP. Rituximab:a review of its use in non Hodgkin's lymphoma and chronic lymphocytic leukaemia. *Drugs* 2003; 63:803-843.
20. Bizzaro, N., et al. EDTA-Dependent Pseudothrombocytopenia: A Clinical and Epidemiological Study of 112 Cases, With 10-Year Follow-Up. *American Journal of Hematology* 1995; 50:103-109.
21. Drachman J G. Inherited thrombocytopenia:when a low platelet count does not mean ITP. *Blood* 2004; 103:390-398.
22. Harris V. K. N and Harris S. Harris Platelet Syndrome—

- Underdiagnosed and unrecognized. *Archives of Pathol Lab Med* 2008; 10:1546.
23. Briggs C, Harrison P, Machin SJ. Continuing developments with the automated platelet count. *Int J Lab Hem* 2007; 29:77-91.
 24. Schoori M, schoori M, Oomes J, van Pelt J. New fluorescent method (PLT-F) on Sysmex XN2000 hematology analyzer achieved higher accuracy in low platelet counting. *Am J Clin Pathol* 2013; 140:495-9.
 25. Stasi R. How to approach thrombocytopenia. *Hematology Am Soc Hematol* 2012; 2012:191-197.
 26. Martel N, Lee J, Wells PS. Risk for heparin induced thrombocytopenia with unfractionated heparin and low molecular weight heparin prophylaxis:a meta analysis. *Blood* 2005; 106:2710-2715.
 27. Arber DA, Orazi A, Hasserjian R, Thiele T, et al. The 2016 World health Organisation classification of myeloid neoplasms and acute leukaemias. *Blood* 2016; 127:2391-2405.
 28. Briggs C, Kunka S, Hart D, Ogumi S, et al. assessment of immature platelet fraction in peripheral thrombocytopenia. *British Journal of Haematology* 2004; 126:93-99.
 29. Dadu T, Sehgal K, Joshi M, Khodaiji S. Evaluation of the immature platelet fraction as an indicator of platelet recovery in dengue patients. *Int J Lab Hematol* 2014; 36:499-504.